


1

**Welcome!**

Attach your Watch Receiver to the charger, and connect it to an USB outlet. Let it charge for at least 15 min.

Remove your watch from the charger, select your language and your watch is ready for pairing.




2

**Pair your watch**

- ◀ First, make sure your bridge is installed correctly. Then, hold down the top test button until the front LEDs blink alternately to initiate pairing.

A screen appears to confirm when your watch has been paired. ▶



3

**Test the connection**

- ◀ Press the bridge top test button to send a signal to your watch.

A status screen confirms that your watch is connected. ▶



4

### Position your watch

Wear your watch comfortably on your wrist, yet firmly enough to allow it to generate distinct vibrations when you receive notifications.



5

### Wake-up



Raise your wrist or press any of the side keys to wake up the screen.

### Navigation



Tap the screen to e.g. open a panel.



### Watch face



Indicates the battery level.



Displays missed notifications.

Tue 14

Shows the day and date.

6

### Home key



Press to go back to the watch face.



Press and hold to turn your watch on or off.



### Back key



Press to open the menu or go back one step.



Press and hold for 3s to call for attention.

### Discover all features!



[bellman.com/watch-guide](https://bellman.com/watch-guide)

4

### Position your watch

Wear your watch comfortably on your wrist, yet firmly enough to allow it to generate distinct vibrations when you receive notifications.



5

### Wake-up



Raise your wrist or press any of the side keys to wake up the screen.

### Navigation



Tap the screen to e.g. open a panel.



### Watch face



Indicates the battery level.



Displays missed notifications.

Tue 14

Shows the day and date.

6

### Home key



Press to go back to the watch face.



Press and hold to turn your watch on or off.



### Back key



Press to open the menu or go back one step.



Press and hold for 3s to call for attention.

### Discover all features!




[bellman.com/watch-guide](https://bellman.com/watch-guide)

1

**Welcome!**

Attach your Watch Receiver to the charger, and connect it to an USB outlet. Let it charge for at least 15 min.

Remove your watch from the charger, select your language and your watch is ready for pairing.




2

**Pair your watch**

- ◀ First, make sure your bridge is installed correctly. Then, hold down the top test button until the front LEDs blink alternately to initiate pairing.

A screen appears to confirm when your watch has been paired. ▶



3

**Test the connection**

- ◀ Press the bridge top test button to send a signal to your watch.

A status screen confirms that your watch is connected. ▶



## Position your watch

Wear your watch comfortably on your wrist, yet firmly enough to allow it to generate distinct vibrations when you receive notifications.



## Wake-up



Raise your wrist or press any of the side keys to wake up the screen.

## Navigation



Tap the screen to e.g. open a panel.



## Watch face



Indicates the battery level.



Displays missed notifications.

Tue 14

Shows the day and date.

## Home key



Press to go back to the watch face.



Press and hold to turn your watch on or off.



## Back key



Press to open the menu or go back one step.



Press and hold for 3s to call for attention.

## Discover all features!



BE3330\_048MAN4.0

[bellman.com/watch-guide](https://bellman.com/watch-guide)

4

### Position your watch

Wear your watch comfortably on your wrist, yet firmly enough to allow it to generate distinct vibrations when you receive notifications.



5

### Wake-up



Raise your wrist or press any of the side keys to wake up the screen.

### Navigation



Tap the screen to e.g. open a panel.



### Watch face



Indicates the battery level.



Displays missed notifications.

Tue 14

Shows the day and date.

6

### Home key



Press to go back to the watch face.



Press and hold to turn your watch on or off.



### Back key



Press to open the menu or go back one step.



Press and hold for 3s to call for attention.

### Discover all features!



[bellman.com/watch-guide](https://bellman.com/watch-guide)