

Hearing and Alerting Solutions for Veterans



Program Overview

Contract Details

2023 - 2028 VA Assistive Listening Devices Contract
In partnership with ADCO Hearing Products



For Sales and General Info

Email: sales@adcohearing.com

Phone: (800) 726-0851

Service and Resources

Useful Resources

Resources for all products, including user manuals, datasheets, and instructional videos, are easily accessible online for both provider and patients.

Training and Continuing Education

From in-person demonstrations to online virtual sessions, we are pleased to offer flexible training and education modules. To schedule your training, contact us or book online. AudiologyOnline CEU courses are also available.

Warranty

A two-year manufacturer's warranty covers main products against defects in materials and workmanship. For full terms and details, refer to the user manuals.

Provider and Patient Customer Support

Provider Email: va.support@bellman.com

Patient Email: support.bsn@bellman.com

Phone: (833) 235-5626

Hours: 8:00am - 8:00pm EST Monday - Friday



programs.bellman.com/us/va

Autonomy through Technology

At Bellman & Symfon, we design for more than just function. We create reliable, user-friendly products that help people with hearing challenges stay safe, independent, and connected.

Founded in 1989, Bellman & Symfon has partnered with hearing care providers, assistive technology professionals, and the wider hearing community to develop innovative solutions that are intuitive, dependable, and life-changing.



Bellman Personal Amplifiers



Maxi Pro Mino

Discover the Joy of Better Hearing

Bellman & Symfon’s Maxi Pro and Mino are fully digital, easy-to-use personal amplifiers designed to make conversations, TV audio, and phone calls clearer and more enjoyable in any environment. Both models work with wired headphones, wired earbuds, or hearing aids*, and use advanced digital signal processing to reduce background noise while preserving natural speech clarity.

*Bellman Neckloop (BE9159) is required to use with hearing aids equipped with T-coil functionality.

Maxi Pro Personal Amplifier

BE2021-P01



Bluetooth enabled, for amplified personal conversations and cell phone calls.

A Bluetooth digital listening solution that helps you stay connected. Amplify personal conversations, take phone calls on your smartphone, receive audio from Bluetooth devices, and stream TV audio [†].



Fully Digital
Clear and noise free sound



Omnidirectional Mic
Sound coverage in all directions



Bluetooth
Connects wirelessly for audio streaming



TV Streaming [†]
Streams TV sound in full stereo

[†] Streaming TV audio to Maxi Pro requires the Bellman TV Streamer (BE2022), available for purchase through Prosthetics.



Scan to learn more

Mino Personal Amplifier

BE2030-P01



Directional microphone allows you to point and capture speech, even in noisy environments.

A digital personal amplifier that helps you hear better in most everyday situations. Features a directional microphone to focus on speech, an omnidirectional microphone for general listening, and a built-in telecoil.



Fully Digital
Clear and noise free sound



Omni + Directional Mics
Handles both quiet and busy surroundings



Telecoil
Compatible with public venues



Slim & Lightweight
Fits discreetly in your pocket



Scan to learn more

Bellman Alarm Clocks



Vibio Digital Alarm Clock

Effortless Wake-Up Solutions

For people who are heavy sleepers or have hearing loss, a traditional alarm clock may not be enough. Bellman & Symfon’s Digital Alarm Clock and Vibio are user-friendly solutions designed to complement any lifestyle, at home or on the go. Features strong vibrations and dependable backup functions during emergencies.

Digital Alarm ClockBE1380-P01



Wakes you up with increasing sound, bright flashing lights, and vibrations.

A digital alarm clock with ascending audio up to 100dB with multi-tone frequencies, vibrant visual cues, and strong vibrations. Easy to use, with large dials and intuitive setup. Battery backup included. 12/24-hour time selection.

- Digital Clock**
User-friendly, great ergonomics
- Loud Sound**
Ascending sound up to 100dB
- Intense Lights**
Wakes you with bright flashing lights
- Strong Vibrations**
Bed Shaker included for tactile alarms



Scan to learn more

Vibio Bluetooth Bed ShakerBE1221-P01



Portable alarm bed shaker for home or travel with phone notifications.

Vibio is a wireless bed shaker that connects to your mobile via Bluetooth with a free app for both iOS and Android. In the Vibio app, you can create multiple alarms, choose vibration power and allow notifications from incoming calls or messages.

- App Controlled**
Customize your alarms with the free Vibio app
- Strong Vibrations**
Customizable strength from mild to strong
- Phone Notifications**
Notifies you of incoming calls and text messages
- Rechargeable**
Lithium-ion batteries with power cable



Scan to learn more

Bellman Alerting System



Stay Safe, Independent, and Connected

The simplest and most affordable way for people with hearing challenges to stay safe and independent at home. Whether it's a doorbell, phone call, baby crying, smoke alarm, or carbon monoxide leakage, the Bellman Alerting System provides reliable notifications through sound, lights, or vibrations. A family of transmitters to detect events, and receivers to alert, that work wirelessly as one system. All items work out of the box, no pairing* or WiFi required.

* Bellman Bluetooth Bridge and Watch Receiver require Bluetooth pairing.

Transmitters – Detect events

Fire and smoke



Smoke Alarm Transmitter
BE1284

Carbon monoxide leak



CO Alarm Transmitter
BE1210

Button press



Push Button Transmitter
BE1240

Doorbell ringing



Door Transmitter
BE1411

Landline phone call



Telephone Transmitter
BE1431

Sound trigger



Baby-Cry Transmitter
BE1491

Receivers – Alert types

Bedside Alerts: Flashing lights, strong vibrations, and loud sound



Alarm Clock Receiver (with Bed Shaker)
BE1380-A

Wearable Alerts: Event-based vibrations



Pager Receiver + Charger + Bed Shaker
BE1230 + BE1260 + BE1271

Tabletop Alerts: Bright flashing lights



Flash Receiver
BE1442

Portable Alerts: Adjustable sound up to 93 dB



Portable Receiver
BE1033

On-the-wrist Alerts: Subtle vibrations, and optional smartphone alerts



Watch Receiver Bluetooth Bridge
BE3330 + BE1521

Smartphone Alerts: Bellman Assistant App, free on iOS and Android



Bluetooth Bridge
BE1521

 Bed Shaker (BE1271) compatible



Scan to learn more

PACKAGE EXAMPLES

Safety Solutions

Receive alerts through sound, lights, or vibrations when smoke or carbon monoxide leakage is detected, even while asleep.

Example includes

- BE1284** Smoke Alarm Transmitter
- BE1210** CO Alarm Transmitter
- BE1442** Flash Receiver
- BE1380-A** Alarm Clock Receiver with Bed Shaker



Caregiver Support

Alert a caregiver with movement or the push of a button. Alerts with a Watch Receiver or on your smartphone.

Example includes

- BE1240** Push Button Transmitter
- BE9026** Contact Mat
- BE1521** Bluetooth Bridge
- BE3330** Watch Receiver



Home Alerting

Alerting for guest and calls. Full home coverage with wearable alerts, bright flashing lights, and bedside notifications.

Example includes

- BE1411** Door Transmitter
- BE1442** Flash Receiver
- BE1380-A** Alarm Clock Receiver with Bed Shaker
- BE1521** Bluetooth Bridge
- BE3330** Watch Receiver



Bellman Accessories

Alerting Accessories

Add vibration



Bed Shaker
BE1271

Compatible with

- BE1442** Flash Receiver
- BE1033** Portable Receiver
- BE1260** Pager Charger

Pressure sensor



Contact Mat
BE9026

Compatible with

- BE1431** Telephone Transmitter
- BE1521** Bellman Bluetooth Bridge
- BE1491** Baby Monitor
- BE1380-A** Alarm Clock Receiver

Door / window open



Magnetic Switch
BE9023

Compatible with

- BE1431** Telephone Transmitter
- BE1521** Bellman Bluetooth Bridge
- BE1491** Baby Monitor
- BE1380-A** Alarm Clock Receiver

Easy smoke / CO installation



Mounting Kit
BE9175

Compatible with

- BE1284** Smoke Alarm Transmitter
- BE1210** CO Alarm Transmitter

Personal Amplifier Accessories



Maxi Pro
Headphone /w Mic
BE9233



Mino Stereo
Headphones
BE9122



Maxi Pro
Earphones /w Mic
BE9229



Mino Stereo
Earphones
BE9214



Bellman
Neckloop
BE9159



Mino
External Mic
BE9409



Maxi Pro
TV Streamer*
BE2022



Mino TV
Cable Kit
BE9126

*TV Streamer (BE2022) available for purchase through Prosthetics.



Our Mission

We are driven to make a difference for people with hearing challenges. We realize this by offering solutions that empower people to stay safe, independent, and socially active.

Bellman & Symfon is proud to offer products under the 2023–2028 VA ALD Contract, in partnership with ADCO Hearing Products.

Bellman & Symfon North America Inc.

Web: programs.bellman.com/us/va

Email: va.support@bellman.com

Phone: (919) 439-8383

Toll-free: (833) 235-5626



programs.bellman.com/us/va