

Driven to Make a Difference

For over 35 years

As a global leader in assistive listening and alerting technology, Bellman & Symfon is proud to offer products in the 2023-2028 VA ALD Contract, in partnership with ADCO Hearing Products.



HEARING AMPLIFIERS

	Part Number	ROES Name	Sound Processing	Connectivity	Microphone	Rechargeable Batteries	Typical Applications
	BE2021-P01	Bellman Maxi Pro Amplifier	Digital	Bluetooth	Omnidirectional	Included, 70hr/charge	• 1-to-1 or small group setting • Cell phone amplification • Stream tablet, computer audio
	BE2030-P01	Bellman Mino Amplifier	Digital	Telecoil	Omnidirectional & Directional	Included, 29hr/charge	• 1-to-1 or small group setting • Listen in noisy environments • Connect to Loop System

ACCESSORIES

Part Number	ROES Name	Typical Application	Part Number	ROES Name	Typical Application
 BE9233	Maxi Pro Headphones /w Mic	Comfortable on-ear headset with mic.	 BE9122	Mino Stereo Headphones	Comfortable stereo on-ear headset.
 BE9229	Maxi Pro Earphones /w Mic	Earbuds with mic, limits sound loss.	 BE9124	Mino Stereo Earphones	Earbuds, limits sound loss.
 BE9159	Bellman Neckloop	Connects amplifiers to hearing aids equipped with T-coil.	 BE9409	Bellman Mino External Mic	External mic for meetings, TV listening, and more.
 BE9147	Maxi Pro/Mino Travel Case	Protective travel case for Maxi Pro and Mino.	 BE9126	Bellman Mino TV Cable Kit	16 feet long 3.5mm cable to connect Mino to TV
 BE2022	Maxi Pro TV Streamer	Available for purchase through VA Prosthetics.			

DIGITAL ALARM CLOCK / WIRELESS BED SHAKER

	Part Number	ROES Name	Bluetooth	Alert Types	Power	Typical Applications
	BE1221-P01	Wireless Alarm Clock/Bed Shkr BT (Vibio)	Yes	Adjustable vibrations, set via free Vibio App	Rechargeable	• Portable bed shaker • Travel alarm • Alerts to incoming calls or texts
	BE1380-P01	Bellman Alarm Clock W/Bed Shaker	No	Sound, flashing light, and/or vibrations	Power cord with backup battery	• Bedside alarm clock



SOLUTIONS BY BELLMAN & SYMFON
PROGRAMS.BELLMAN.COM/US/VA
VA.SUPPORT@BELLMAN.COM | 833.235.5626

ALERTING SYSTEMS

- A customizable system of transmitters and receivers that connects as many devices as needed.
- A modular system, ready to use out of the box with no Wi-Fi required.

TRANSMITTERS - DETECT EVENTS

	Part Number	ROES Name	Battery Life (Up to)	Typical Applications
	BE1284	Bellman Smoke Alarm Transmitter, UL certified	10 Years	Detects fire or smoke at an early stage.
	BE1210	Bellman Carbon Monoxide Alarm Transmitter, UL certified	5 Years	Detects CO leakage at different levels.
	BE1240	Bellman Push Button Transmitter	2 Years	Works as a push-button doorbell. Also works as a call-for-attention button.
	BE1411	Bellman Door Transmitter	5 Years	Detects the sound of an existing door chime.
	BE1431	Bellman Telephone Transmitter	5 Years	Detects landline phone calls. Also works as a multi-purpose connector for accessories below.
	BE1491	Bellman Baby-Cry Transmitter (Sound Detector)	6 Months	Detects baby crying. Also works as a sound detector for other sounds.

RECEIVERS - ALERT TYPES

	Part Number	ROES Name	Placement	Alert Type			Typical Application
				Sound	Lights	Vibration	
	BE1380-A	Bellman Receiver w/ Bed Shaker, UL certified (Alarm Clock Receiver)	Bedside				Bedside alerts with sound, flashing lights, and strong vibrations . Bed Shaker and backup battery included.
	BE1230	Bellman Pager Receiver	Wearable				Wearable alerts with event-based vibration patterns , and LED event display panel. Bed Shaker compatible, with Pager Charger.
	BE1442	Bellman Flash Receiver, UL certified	Tabletop				Alerts with bright flashing lights . Rotating top directs flashes. Backup battery included. Bed Shaker compatible.
	BE1033	Bellman Portable Receiver	Tabletop				Portable alerts with adjustable sound , up to 93 dB. Bed Shaker compatible.
	BE3330	Bellman Bluetooth Watch Receiver	Wearable				Provides discreet vibrations on the wrist . Requires the Bellman Bluetooth Bridge.
	BE1521	Bellman Bluetooth Bridge	Mounted				Bluetooth Transceiver. Enables alerts on cell phone via the free Bellman Assistant app; enables connection with Watch Receiver .

ACCESSORIES - EXPAND

Part Number	ROES Name	Typical Application	Part Number	ROES Name	Typical Application
	BE1260	Bellman Pager Charger		BE9026	Bellman Contact Mat
		Table-stand charger for BE1230 Pager Receiver.			Triggers an alert when someone steps on the mat.
	BE1271	Bellman Bed Shaker, UL certified		BE9023	Bellman Magnetic Switch
		Compatible with BE1380-A, BE1442, BE1033, BE1260.			Triggers an alert when a door or window opens.
	BE9175	Bellman Mounting Kit		BE9250	Bellman Mobile Phone Sensor
		Quick and easy CO & Smoke Alarm mount.			Detects cell phone and tablet notifications.



SOLUTIONS BY BELLMAN & SYMFON
PROGRAMS.BELLMAN.COM/US/VA
VA.SUPPORT@BELLMAN.COM | 833.235.5626

